

Ordering Guide

To place your order call: 1-877-665-4381



- Patient Height _____ (Less than 6.0" recommend Independence II / Greater than 6.0" recommend Independence IV)
- Patient Weight _____ (Under 250 pounds recommend Independence II / Greater than 250 pounds recommend Independence IV)
- Does the patient plan to travel with the toilet lift? Y | N (If yes, consider locking caster accessory)
- Does the patient have an outlet in the bathroom? Y | N (If no, consider battery kit option)
- Does patient live in an area prone to power outages? Y | N (If yes, we recommend battery back-up kit)
- Will the product be moved between the bedroom and the bedside frequently? Y | N (If yes, consider locking casters and commode bucket)
- Does the patient's medical condition limit the use of their fingers? Y | N (If yes, consider big button hand pendant and/or foot switch)
- Would the patient prefer a soft padded toilet seat? Y | N (If yes, we recommend a soft padded toilet seat)
- Does the patient have difficulty performing the hygienic cleaning after toileting? Y | N (If yes, we recommend the electric washing bidet)

Item Number	Item Description	Item Description 2	Price	TOTALS
<input type="checkbox"/> LS300-X-101	LifSeat Independence II	Independence II unit AC powered with removeable washable splash guard	\$1,299.00	\$
<input type="checkbox"/> LS400-X-100	LiftSeat Independence IV	Independence IV unit AC Powered with removeable washable splash guard	\$1,599.00	\$
<input type="checkbox"/> LS-BRN-1000	Brondell Swash 1000	Electric Toilet Seat Washing Bidet	\$649.00	\$
Optional Accessories				
<input type="checkbox"/> LS400-I-EXT	Leveling Feet, Extended	For use with Independence IV - Set of 4 Height Extenders	\$74.00	\$
<input type="checkbox"/> LS900-TS-100	Soft Padded Toilet Seat	For use with Independence II or IV models	\$45.00	\$
<input type="checkbox"/> LS000-CB-100	Reusable Drop-in Commode Bucket	Converts over the toilet into bedside commode	\$24.00	\$
<input type="checkbox"/> LS300-HH-500-I2	Big Button Hand Pendant	For use with Independence II	\$135.00	\$
<input type="checkbox"/> LS300-HH-500-I4	Big Button Hand Pendant	For use with Independence IV	\$135.00	\$
<input type="checkbox"/> LS300-J90035	Foot Switch	For use with Independence II and IV	\$205.00	\$
<input type="checkbox"/> LS900-BK-101	Battery Kit	Includes: Battery, Battery Charger, Mounting Bracket w/Extension Cord	\$300.00	\$
<input type="checkbox"/> LS900-I2W-100	I2 Locking Medical Grade Caster	Set of 4 Medical Grade Caster for Independence II	\$90.48	\$
<input type="checkbox"/> LS900-I4W-100	I4 Locking Medical Grade Caster	Set of 4 Medical Grade Caster for Independence IV	\$90.48	\$
<input type="checkbox"/> LS900-INST-HOME	Home Assembly and Installation	Assembly and Installation of Toilet Lift - Pricing varies on distance from major city and market labor costs - estimate \$250.00	TBD	
<input type="checkbox"/> LS910-INST-HOME BIDET	Home Assembly and Installation	Assembly and Installation of Electric Washing Bidet to toilet (does not include installation of an electrical outlet by the toilet)	TBD	
<input type="checkbox"/> I2 Shipping	Shipping Cost of Independence II	Shipping Cost - if product is returned, customer will be refunded full price including shipping and LiftSeat will pay the return shipping cost.	\$45.00	\$
<input type="checkbox"/> I4 Shipping	Shipping Cost of Independence IV	Shipping cost for the Independence IV model - if product is returned, customer will be refunded full price including shipping and LiftSeat will pay the return shipping cost.	\$65.00	\$
NOTES:			TOTAL	\$

Warranty: 2 Years (Metal Frame, Toilet Seat) 1 Year (Paint, Electric Control Box, Hand Pendant, Motor, Battery Charger, Battery)



“LiftSeat is really a key to increased independence and increased safety for many of our patients”

C. Neumann, OT Manager

To place your order call: 1-877-665-4381

COMMON CONFIGURATIONS



**Independence II
with electric wash bidet**



**Independence II
with locking wheels**



**Independence IV
with electric wash bidet
and raised seating height extenders**

ACCESSORIES



Locking Wheels



**Raised Seating
Height Extenders**



Hand Pendant Options



Foot Pedal Option



Back up battery



Soft Toilet Seat



Drop In Bucket
For Bedside Commode Use



Electric Washing Bidet

Other configurable options include battery only powered units, splash guards and a back rest for bedside commode use

Assembly and Installation

1. Screw in leveling feet, attach toilet seat and insert handle bars
2. Remove existing toilet seat from toilet
3. Align LiftSeat over toilet bowl opening
4. Plug into GFI Outlet
5. Test range of motion (down / up)

OR

LiftSeat conducts white glove service of assembly and installation = \$TBD

“Recently purchased one for my home use. It has given me back independence and privacy when using toilet. The company was great to work with.

I highly recommend”

B. Rosenberg, LiftSeat user